



People.Health.Care.

Sentina® Ivory Latex Powder-free examination glove



Sentina® Ivory Latex Powder-free examination glove



Produktinformationen:

- High wearing comfort and excellent tactile sensitivity thanks to optimised grammage.
- Secure grip thanks to textured finger tip.
- With synthetic inner coating for quick and easy donning.

Fields of application:

- For use in hospitals and in outpatient settings.
- For professional cleaning.
- Suitable for laboratories, pharmaceutical companies and for contact with food.

Sentina® Ivory Latex

Powder-free examination glove,
non-sterile, off-white

| Size | REF | Shipping units (pieces) |
|------|---------|-------------------------|
| XS | 163 040 | 100/1000 |
| S | 163 041 | 100/1000 |
| M | 163 042 | 100/1000 |
| L | 163 043 | 100/1000 |
| XL | 163 044 | 90/900 |

Product features:

| | |
|---|--|
| Product type: | Powder-free examination and disposable protective gloves, non-sterile |
| Material: | Natural rubber latex |
| Colour: | Off-white |
| Glove shape: | Flat with rolled rim, ambidextrous |
| Glove surface: | Textured finger tip |
| Size/ Overall length: | XS, S, M, L, XL / median 240 mm (according to EN 455-2) |
| Wall Thickness: (single wall thickness) | Palm area: typ. 0.10 mm Finger: typ. 0.12 mm |
| Barrier properties: | AQL 1.5 (according to EN 455-1) |
| Tear resistance: | Median \geq 6 Newton (according to EN 455-2) |
| Shelf life: | 5 years (in original packaging when stored in accordance with ISO 2230) |
| Packaging: | 100/90 pcs per dispenser box, 10 dispenser boxes per transport carton |
| Labeling and conformity to standards: | EN 420, EN ISO 374-1, EN ISO 374-2, EN ISO 374-3, EN 374-4, EN ISO 374-5, EN 421, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D3578, EN 16523-1 |
| Intended purpose: | Disposable medical gloves class I according to Medical Device Regulation (EU) 2017/745. Disposable protective gloves category III (temporary protection against chemical effects)* according to PPE regulation (EU) 2016/425. Suitable for contact with food according to regulation (EG) 1935/2004. |
| Additional information: | This product contains natural rubber latex and chemical accelerators. In case of intolerance, do not use these gloves. A list of the substances contained in the glove is available on request. |

* Please note that product must be used in accordance with its chemical resistance and the respective declaration of conformity. Failure to observe this information, in particular with regard to chemical resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Lohmann & Rauscher is not liable for incorrect use of the gloves nor will Lohmann & Rauscher be liable for improper storage and handling. In case of indications of damage or unusual signs of wear and tear to the gloves, the gloves should not be used. In case of any doubt, expert advice should be sought before use. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments.